

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS22 4730

## PsaD | PSI-D subunit of photosystem I (algal)

## **Product information**

Immunogen PsaD of Chlamydomonas reinhardtii, UniProt: Q39615

**Host** Rabbit

Clonality Polyclonal

Purity Serum

Format Lyophilized

Quantity 50 μl

**Reconstitution** For reconstitution add 50 μl of sterile water

Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles, Please, remember to spin tubes briefly prior to opening them to avoid any losses that might occur from lyophilized

material adhering to the cap or sides of the tubes

## **Application information**

**Recommended dilution** 1: 5000 - 1 : 10 000 (WB)

Expected | apparent 21.3 kDa

Confirmed reactivity Chlamydomonas reinhardtii

**Additional information** Can be provided with ProClin, if requested