

This product is **for research use only** (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91112 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no **AS06 109**

PsaN | PSI-N subunit of photosystem I

Product information

Immunogen	KLH-conjugated synthetic peptide chosen from <i>Arabidopsis thaliana</i> PsaN protein sequence At5g64040
Host	Rabbit
Clonality	Polyclonal
Purity	Immunogen affinity purified serum in PBS pH 7.4.
Format	Lyophilized
Quantity	200 µg
Reconstitution	For reconstitution add 133 µl of sterile water
Storage	Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.

Application information

Recommended dilution	1 : 1000 (WB)
Expected apparent MW	9,7 kDa
Confirmed reactivity	<i>Arabidopsis thaliana, Hordeum vulgare, Nicotiana tabacum</i>
Predicted reactivity	<i>Catalpa bungei, Phaseolus vulgaris, Oryza sativa, Populus balsamifera, Zea mays, Vitis vinifera</i> Species of your interest not listed? Contact us
Not reactive in	<i>Chlamydomonas reinhardtii, Synechococcus</i> sp. PCC7942
Selected references	Hansson et al. (2007). Knock-out of the chloroplast encoded PSI-J Subunit of Photosystem I in Nicotiana tabacum: PSI-J is required for efficient electron transfer and stable accumulation of photosystem I. FEBS J. 274: 1734-1746.